

Welcome to the world of real-time CGM

Jennifer Ruck, BSN, RN, CDCES Sr. Clinical Account Manager jennifer.ruck@Dexcom.com





LBL016304 Rev00


*If your glucose alerts and readings from the G6 do not match symptoms or expectations, use a blood glucose meter to make diabetes treatment decisions.

Agenda


- What is CGM?
- G6 components
- G6 app and Follow app
- Sensor insertion
- Treatment decisions
- Troubleshooting
- Tech Support

LBL016346 Rev 001

BGM vs. CGM



BGM vs. CGM



Real-Time CGM is the standard of care for patients on intensive insulin therapy.





Recognized as a standard of care in diabetes management by the American Diabetes Association (ADA), American Association of Clinical Endocrinologists (AACE) and the Endocrine Society.^{1,2,3}

Continuous glucose monitoring (CGM) use has proven to both reduce A1C and time spent in hypoglycemia regardless of delivery method.^{4,5}

Regular CGM use is also associated with an increase of time spent in range,^{4,6} a reduction of severe hypoglycemic incidents,^{7,8} and associated with many behavioral changes that promote significant positive change in diabetes self-management.⁹

1 American Diabetes Association. (2016). Glycensic Targets. Standards of Medical Care. Diabetes Care, 533-540

¹Formace Y, Grunberger G, Anhale H et al. COMPRINGUES GUICOSE MONITORING: A CONSERVENCE OF THE AMERICAN ASSOCIATION OF CLINICAL INDOCENNOLOGISTS AND AMERICAN COLLEGE DE ENDOCENNOLOGY. Endoor Pract. 2009;22 [j1:008-1021. ³Peters A, Almann A, Battelino T et al. Diabetes Technology – Continuous Subcutaneous Imulin Indusion Theory and Continuous Subcutaneous Imulin Adults: An Endoorine Society Clinical Practice Guidelines. J Clinical Endoorino (Metab. 2006;20:06-2034. ³Sougal.), Petrubilizavia., Head M et al. Comparison of Phirent Treatment Modulitation Trype 1. Diabetes, Including Sensor-Augmented Intuin Regiment, in 52 Weeks of Follow-Up: A CONISAR Study. Diabetes Technology & Theorpeutics. 2016;21(0):532-538.

*Batteling T.P.M., Brating N., Nimri R., Oskansson P., Bolinder J. Effect of Continuous Glucose Monitoring on Hypoglycersia in Type 1 Diabetes: Diabetes: Care. 2011;54(4):725-800.

⁴Beck RW, Hiddlesworth T, Ruedy K, et al. Effect of continuous glucose monitoring on glycemic control in adults with hype 1 diabetes using insulin injections: The diamond randomized clinicalized. JAMA. 2017;317(4):171-378. doi:10.1001/jama.2016.19375. ¹Choudhary P, Ramanamy S, Green L et al. Real-Time Continuous Glucose Monitoring Sign Boardly Reduces Severe Hypoglycemia in Hypoglycemia-Unasere Patients With Type 1 Diabetes. Diabetes Care. 2013;30(2):100-4162.

^aLind M, Pdomby, W, Hinch, L, et al. Continuous Glaccos Monitoring in Convestional Therapy for Gloomic Control in Adults With Type 1 Diabetes Treated With Multiple Daily Injections – The GDLD Randomized DivisitTrial. [published January 2017]. JAMA

⁰Haas L, Marjin Lik M, Beck J, et al. National Standards for Diabetes Self-Management Education and Support, Diabetes Care, Volume 35, November 2012

APPROVED FOR 2+ YEARS

The only CGM FDA-approved for ages two years and older

Holden and Zola, T1Ds, Dexcom Warriors







* For a list of compatible devices see: dexcom.com/compatibility

** In this guide, we will only be discussing the smart device option. For receiver information see the Dexcom G6 User Guide (dexcom.com/guides).



Dexcom G6 App

Dexcom Clarity App

*For a list of compatible devices see: dexcom.com/compatibility


SMART DEVICE USE

If the student uses a smart device it will need to be with them at all times and the volume will need to be loud enough to hear low and high alerts.

The 504 plan will outline the use of the smart device in a school setting.



LBL015708 Rev 001



TREND ARROWS

Trend arrows shows a student's glucose speed and direction.



DEXCOM FOLLOW®

To receive the student's glucose information, Followers will need to install the Dexcom Follow app on a compatible smart device* and have internet connectivity. When someone is added as a Follower they will be sent an email with instructions.

Any diabetes treatment decision should be based on the student's G6, not information from the Dexcom Follow app.



Dexcom Follow App

* For a list of compatible devices see: dexcom.com/compatibility

TREATMENT DECISIONS

A student may come to you to treat a low or high glucose level. The G6 can be used to make treatment decisions such as eating for a low glucose level or dosing insulin for a high glucose level. Base treatment decisions on the student's 504 plan. **201** mg/dL

Use a meter:


- · If a number and an arrow are not shown on the G6 app
- Any time the student's symptoms or expectations don't match readings. For example, if a student says they feel low but the G6 shows them in a normal range.




LOW AND HIGH ALERTS

Each student should have a set low and high glucose alert.

The display device will either vibrate or beep based on the student's alert settings.





RESPONDING TO ALERTS AND ALARMS

Based on the 504 plan a student may come to you to help respond to alerts and alarms.

The steps you should take are:

- Go into the G6 app on the student's smart device
- Tap OK to clear the alert
- Take action based on the glucose information shown in the G6 app and 504 plan

SENSOR INSERTION


You may need to help a student with sensor insertion. Sensor insertion should be outlined in the 504 plan.

Find sensor insertion video here: dexcom.com/training-videos



TROUBLESHOOTING

A student may come in for help with a screen they have never seen before. Use a meter for treatment decisions if an error is shown.



https://www.dexcom.com/contact



Global Technical Support

Product troubleshooting or replacement inquiries **1-844-607-8398** Available 24 hours a day; 7 days a week

Request a call back from a representative Submit a Patient Support Request Request Sensor Overpatches Chat Live with Dexcom Tech Support (At this time, live chat is only available to US customers, Monday-Sunday from 5AM-5PM PST.)

EVERY STUDENT IS UNIQUE

If other issues come up regarding the student and their CGM, speak to the main caregiver and adjust the 504 plan as needed.

